



YOU CAN VIEW IT ANYTIME, ANYWHERE*



Patients can easily share their glucose data
for more productive diabetes management
conversations with you remotely


FreeStyle
Libre
FLASH GLUCOSE MONITORING SYSTEM


life. to the fullest.®
Abbott

The FreeStyle Libre flash glucose monitoring system is indicated for measuring interstitial fluid glucose levels in adults aged 18 years and older with diabetes mellitus. Always read and follow the label/insert. The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app.
* The LibreView website is only compatible with certain operating systems and browsers. Please check www.LibreView.com for additional information.

Patients can check their glucose with the FreeStyle LibreLink app,* and you see their data with LibreView†




Convenient for your patients

Patients can easily check their glucose using the FreeStyle LibreLink app, their FreeStyle Libre reader, or both




LibreView



Convenient for you

Patients' glucose data are automatically uploaded‡ from their phone to LibreView – no need for office staff to download the data



* The FreeStyle LibreLink app and the FreeStyle Libre reader have similar but not identical features. A finger prick test using a blood glucose meter is required during times of rapidly changing glucose levels when interstitial fluid glucose levels may not accurately reflect blood glucose levels or if hypoglycemia or impending hypoglycemia is reported by the FreeStyle LibreLink app or when symptoms do not match the app readings. The FreeStyle Libre sensor communicates with the FreeStyle Libre reader that started it or the FreeStyle LibreLink app that started it. A sensor started by the FreeStyle Libre reader will also communicate with the FreeStyle LibreLink app. The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink requires registration with LibreView.

† The LibreView website is only compatible with certain operating systems and browsers. Please check www.LibreView.com for additional information.

‡ Use of FreeStyle LibreLink requires registration with LibreView. Automatic upload requires a wireless internet connection or mobile data connection.

Generate glucose reports on demand to support effective diabetes management



Easy upload

Glucose data are automatically uploaded* from patients' FreeStyle LibreLink app to LibreView and can be manually uploaded from the FreeStyle Libre reader



Remote access

Log in and view standardized reports from any internet-connected computer



Secure and reliable

Glucose data are available online through a secure, cloud-based system at no cost




Note: Select LibreView reports shown here.
Complete reports available at LibreView.com.

Patients can use the FreeStyle LibreLink app to easily share their glucose data with you remotely




CONNECT TO A PRACTICE using the FreeStyle LibreLink app

Connect to a Practice is a FreeStyle LibreLink feature that lets patients easily share their glucose data with your practice through LibreView



A secure online tool to view patients' glucose data

Your patients' glucose data are automatically uploaded* to LibreView, a cloud-based system that enables you to access reports on demand



LibreView



Easily connect

Patients can initiate the connection or they can accept an invitation from your practice

Convenient for you


Once connected, patients' glucose data are automatically uploaded* from the FreeStyle LibreLink app to LibreView – no need for office staff to download the data



Virtual office visits

Easy access to glucose reports for more informed treatment decisions and productive conversations

Easy to invite patients to automatically share their glucose data with your practice

Option 1: Invite patients to connect from LibreView

1. Search for a patient or click the “Add Patient” icon  at the top of the navigator bar
2. In the patient’s profile, select invite and use the email associated with their account
3. Patients receive invitations in their FreeStyle LibreLink app and in their email



Patient accepts invitation in the FreeStyle LibreLink app or in their email





Need help? Call Abbott’s HCP Customer Care line at 1-844-610-1001 or visit the “Support” section on www.LibreView.com

Give patients your LibreView Practice ID so they can easily enable data sharing from the FreeStyle LibreLink app

Option 2: Patients enable data sharing with your practice

Provide your patients with your Practice ID, which is available in LibreView under “My Practices”



Patient enters the Practice ID in the FreeStyle LibreLink app to enable connection




Need help? Call Abbott's HCP Customer Care line at 1-844-610-1001 or visit the “Support” section on www.LibreView.com



YOU CAN VIEW IT

ANYTIME, ANYWHERE*



Patients can share their glucose data for more productive conversations with you remotely

- It's easy for patients to share their glucose data with your practice through LibreView at no cost
- Glucose data are automatically uploaded† from patients' FreeStyle LibreLink app to LibreView and can be manually uploaded from the FreeStyle Libre reader
- Use LibreView reports to see your patients' glucose patterns and trends to support effective diabetes management



**Recommend the FreeStyle LibreLink app
and easily connect with your patients through LibreView**



* The LibreView website is only compatible with certain operating systems and browsers. Please check www.LibreView.com for additional information.

† Use of the FreeStyle LibreLink app requires registration with LibreView. Automatic upload requires a wireless internet connection or mobile data connection.

© 2021 Abbott. FreeStyle, Libre, and related brand marks are marks of Abbott. Product images are for illustrative purposes only. Other trademarks are the property of their respective owners.

ADC-21075 v2.0