

# FreeStyle Libre Brand Compatibility Guide



FreeStyle Libre system		FreeStyle Libre 2 system
<b>SENSOR</b>		
FreeStyle Libre sensor PIN: 97799171		FreeStyle Libre 2 sensor PIN: 97799075
<b>COMPATIBLE MOBILE APP</b>		
FreeStyle LibreLink		FreeStyle Libre 2 app
<b>COMPATIBLE READER</b>		
FreeStyle Libre reader PIN: 97799170		FreeStyle Libre 2 reader PIN: 97799074
<b>INDICATED AGES</b>		
18 years and older		4 years and older
<b>GLUCOSE ALARMS</b>		
None available		Three optional alarms – low glucose, high glucose & signal loss
<b>COVERAGE &amp; DISPENSING</b>		
Covered by most private insurers.*		Check your patients' insurance coverage. FreeStyle Libre 2 and FreeStyle Libre prescriptions are not interchangeable.
<b>RETURNS &amp; EXCHANGES</b>		
No store returns or exchanges. Direct customer inquiries to the Abbott Customer Care line: 1-888-205-8296		No store returns or exchanges. Direct customer inquiries to the Abbott Customer Care line: 1-888-205-8296
<b>WHOLESALE ORDER CODES (SENSOR)</b>		
McKesson: TBD K&F: TBD Matrix: TBD PJC: TBD		McKesson: TBD K&F: TBD Matrix: TBD PJC: TBD
<b>WHOLESALE ORDER CODES (READER)</b>		
McKesson: TBD K&F: TBD Matrix: TBD PJC: TBD		McKesson: TBD K&F: TBD Matrix: TBD PJC: TBD



The FreeStyle Libre 2 flash glucose monitoring system is indicated for measuring interstitial fluid glucose levels in people aged 4 years and older with diabetes mellitus. The FreeStyle Libre flash glucose monitoring system is indicated for measuring interstitial fluid glucose levels in adults aged 18 years and older with diabetes mellitus. Always read and follow the label/insert.

The FreeStyle Libre 2 app and the FreeStyle Libre 2 reader have similar but not identical features. Finger pricks are required if readings do not match symptoms or expectations. The FreeStyle Libre 2 sensor communicates with the FreeStyle Libre 2 reader that started it or the FreeStyle Libre 2 app that started it. The FreeStyle Libre 2 app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle Libre 2 requires registration with LibreView.

\* Individual private drug plans can vary. Please check with your plan administrator and/or your insurance company.



# MORE PATIENTS CAN DO IT


WITHOUT LANCETS\*

When dispensing the FreeStyle Libre 2 system,  
here's what to look for:

FreeStyle Libre 2  
product name



FreeStyle Libre 2 reader  
PIN: 97799074



FreeStyle Libre 2 sensor  
PIN: 97799075



The FreeStyle Libre 2 system and FreeStyle Libre system are not cross-compatible.  
The FreeStyle Libre 2 sensor is only compatible with the FreeStyle Libre 2 app and reader.

Learn more at **EasytoFreeStyle.ca**



\* Scanning the sensor does not require lancets.

© 2020 Abbott. FreeStyle, Libre, and related brand marks are marks of Abbott. Product images are for illustrative purposes only. Other trademarks are the property of their respective owners.

ADC-26016